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### Patient Information Leaflet Revamil

<b>The product</b>
Revamil hydrophilic wound gel
5 g/ 18 g gel for the treatment of wounds with honey

**Product description**

Revamil is a wound gel that consists of pure honey, which has been produced under controlled conditions. The carefully controlled production process guarantees that Revamil does not contain any traces of pesticides and that the same quality can be

reproduced consistently. The honey used to make Revamil contains a high proportion of honey enzymes such as glucose oxidase. The highly concentrated sugars cause some of the wound exudate to mix with Revamil gel. The resulting layer helps to strengthen the antibacterial protection provided by the dressing.

<b>Instructions for use</b>
<ol style="list-style-type: none"><li>Before applying Revamil gel first clean the wound.</li> <li>Spread a thin layer of Revamil on the secondary wound dressing.</li> <li>Revamil can also be applied directly to the wound and then covered with a secondary dressing.</li> <li>Attach the secondary dressing in the normal way.</li> <li>Depending on the nature of the wound, replace the dressing every day or every two days.</li></ol>

**Indications**
Revamil wound gel can be used to treat chronic sores (bed sores and different kinds of ulcers), infected wounds, acute wounds, surgical oncological wounds, radiation induced oncological wounds and minor (first- and second-degree) burns.

**Contra indications**

- Do not use the product if the patient is known to be allergic or hypersensitive to honey. An allergic reaction or hypersensitivity can cause redness, irritation and possibly minor pain in the skin or wound.

**Side effects**

- During the first minutes after application, the product can cause a light stinging sensation. The stinging sensation is probably caused by the low pH and high osmolarity of the honey and has no effect on the performance of the product.
- In the event of unexpected side effects after use, contact your doctor or the manufacturer. Side effects can cause inconvenience to the patient.

**Safety precautions**

- In case the honey is crystalized, warm the (closed) tube in lukewarm water to facilitate release of the honey from the tube. Take care that the tube is closed well, when immersed into water. Honey crystallization is a natural process and does not affect the safety or the performance of the product.
- Do not use if honey is leaking from the inner packaging. Leakage can lead to contamination of the product with bacteria and to wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not use if the tube is damaged as this can lead to contamination of the product with bacteria and to wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Always replace the top on the tube after use. An open tube could lead to contamination with bacteria and wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not allow the tube to come into direct contact with the surface of the wound. Direct contact between the wound/ skin and the uncapped tube could lead to contamination with bacteria and wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Before applying Revamil, dampen the wound with 0,9% NaCl. Application of Revamil Gel on a dry wound could lead to lower performance of the product and slower wound healing.
- Do not apply this product for more than 30 days continuously. Treatment for more than 30 days continuously could lead to inconvenience to the patient, such as irritation and/or minor pain.
- The contents of the tube must be used within 3 months of opening. Longer use could lead to lower functionality and contamination of the product. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not treat different patients with the same tube as this can lead to cross-contamination with bacteria and to wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not use the product after the expiration date as this can lead to lower performance of the product and to inconvenience for the patient.

<b>Storing conditions</b>
Store between 10 and 30°C

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