Wound Management Smith & Nephew Medical Ltd. 101 Hessle Road Hull HU3 2BN England T 44 (0) 1482 225 181 F 44 (0) 1482 328 326 www.smith-nephew.com



# Product Safety Data Sheet Bactigras Dressing

## **Approvals**

Author Name: Steven Martin

Signature Date/Time: 12 JAN 2012 11:25 GMTUK Global Date/Time: 12 JAN 2012 11:25 UTC

Representation: Author Signature

Signature Reason: Document Management: Status Change

Signatory Name: Caroline Watson

Signature Date/Time: 24 APR 2012 12:31 GMTUK
Global Date/Time: 24 APR 2012 11:31 UTC
Representation: Engineering Signature

Signature Reason: Document Management: Status Change

1. PRODUCT NAME(S):	BACTIGRAS Dressing.

- 2. INTENDED USES: Antiseptic wound contact layer.
- **3. COMPOSITION:** Contains 0.5% w/w Chlorhexidine diacetate in white soft paraffin.

#### 4. PHYSICAL AND CHEMICAL PROPERTIES:

Releases Chlorhexidine acetate. Combustible solid.

#### **5. HAZARD INFORMATION:**

# **Health effects:**

The product contains substances which are:

- a) Irritants to the eye
- b) Irritants to the mucous membranes

Potential allergy to Chlorhexidine.

6. FIRST AID:	
a) Inhalation	Not applicable.
b) Contact with skin	Not applicable.
c) Contact with eyes	Irrigate well with clean warm water or saline solution. Seek medical aid.
d) Ingestion	Do not induce vomiting. Seek medical aid.

#### 7. FIRE AND EMERGENCY MEASURES:

The product is a combustible solid which gives off toxic fumes, of carbon monoxide, carbon dioxide, hydrogen chloride gas and nitrogen oxides when ignited. In case of fire any standard fire extinguisher may be used except water based. Wear self-contained breathing apparatus.

#### **8. ACCIDENTAL RELEASE MEASURES:**

Take up mechanically and dispose of in accordance with section 14.

## 9. HANDLING AND STORAGE PRECAUTIONS:

Store in a cool, dry place (<25°C). Keep away from sources of ignition.

## 10. EXPOSURE CONTROLS/PERSONAL PROTECTION:

Other than areas of treatment, gloves may be worn to prevent frequent or prolonged exposure.

#### 11. STABILITY AND REACTIVITY:

Not applicable.

#### 12. TOXICOLOGICAL INFORMATION::

See information in section 5.

#### 13. ECOLOGICAL INFORMATION:

The product contains a substance toxic to aquatic organisms and may cause long term adverse effects in an aquatic environment. Do not discharge into waste water or waterways.

## 14. DISPOSAL:

Controlled incineration/ landfill; according to local environmental health guidelines. Do not dispose of into waste water or waterways.

**15. TRANSPORT PRECAUTIONS:** Not applicable.

#### **16. REGULATORY INFORMATION:**

Product contains substances with the following index letter labeling:

Xi - Irritant.

N - Dangerous for the environment.

Document Number: 1002495 Document Part: 000
Document Version: 03 Document Status: Released

Page 3 of 6

Product contains substances with the following R-phrase labeling:

R 36 - Irritant to eyes

R 51 - Toxic to aquatic organisms

R 53 - May cause long term adverse effects to the aquatic environment

#### **17. ADDITIONAL INFORMATION:**

Legal category "P" for product manufactured prior to November 1996, product License number 1613/0001. Product manufactured since this date has been sold as a CE marked medical device.

## 18. NAME, ADDRESS AND TELEPHONE NUMBER OF SUPPLIER:

Compliance & Regulatory Department Smith & Nephew Healthcare Healthcare House 101 Hessle Road Hull

HU3 2BN

Telephone: 01482 222200 Fax: 01482 222211

Email Address: advice@smith-nephew.com

19a) REFERENCE NUMBER	PSDS 008.
19b) DATE OF ISSUE	January 2012.

This information is provided in accordance with the requirements of the UK Health and Safety at Work Act 1974, and specifically in order to assist users of the product to make their 'assessment of health risks' as required by the UK Control of Substances Hazardous to Health Regulation 2002 (COSHH assessments). Provision of this information does not preclude users from seeking advice from other sources as indicated in the COSHH guides. The information is intended to cover potential hazards at the place of work and does not detail medical uses, indications, contra-indications and precautions for the treatment of patient.

## **REASON FOR CHANGE**

Document Version	Section/Paragraph Changed	Change Made	Date
1	Whole document	Migrated from QSI to DMS PDSD 008 rev 5 → DMS Doc. 1002495	25/07/2006
2	NA	Re-format main document and minor amendments COSHH regulations 2002	21/01/2009
3	Whole document	Periodic review: Reformatted to DMS Doc 1003205, no changes made.	12/01/2012

# **RE-VERIFICATION OF DOCUMENT SUITABILITY**

The normal expiry date/time for a printed <u>non-controlled</u> document is 06:00 on the day after the date the document was printed. If such a document is required for more than one day, rather than print a new copy each day, the existing printout can be manually compared to the electronic copy. If it is the same version, you may sign and date the following agreement for each additional day:

"I have compared this paper copy to the original electronic copy on the date shown below and certify that both are the same version number"

Print Name	Signature	Date

Plant Number(s):
None
Deadwat Cada(a)
Product Code(s):
66003650, 66003661, 66007505, 66063650, 66067456, 66067457, 66067461, 66267456
66267457, 66267461, 66333661, 66337457, 66997457, 66997461, 66997505, 7449
7456, 7456G, 7457, 7457G, 7457SGH, 7461, 7461G

\*\*\* END OF DOCUMENT \*\*\*

Document Number:
Document Version:

1002495 03 Document Part: Document Status: 000 Released