

Meliseptol® HBV

Servietter til overflade desinfektion

Indikation

- Overflade desinfektion af alkohol resistent overflader ved medicinske produkter

Egenskaber

- Hurtig og effektiv
- Effektivt på 1 min.
- Mikrobiologisk effektivitet – se rubrik til højre
- Klar til brug servietter
- Fri for aldehyd og alkylaminer
- Dermatologisk testet
- 100 servietter/dispenser
- Mulighed for refill pakke

Produktionbeskrivelse

- Servietten er vædet med Propan-1-ol og Didecyldimethylammonium chlorid
- Servietten er fremstillet af cellulose og syntetiske fibre
- Indeholder parfume



Mikrobiologisk effektivitet

Bacteria (incl. TB bacilli) and yeasts
Fungi (A. brasiliensis)
Enveloped viruses (limited virucidal activity, RKI 1/2004*):
 Incl. Bird flu, HBV, HCV, HIV
Non-enveloped viruses (virucidal acc. RKI 1/2004*):
 Adeno-, Polyoma-, Polio-, Noro-(MNV, FCV) and Rota- viruses
VAH-List (formerly DGHM-List) 5 minutes
 Meliseptol® HBV Tissues are on the List for virucidal products (IHO)

Virketid

Micro-organism	Test Norms	Exposure Time
Surface disinfection Bacteria, yeasts (<i>C. albicans</i>) and enveloped viruses (incl. HBV, HCV, HIV, Vaccinia virus, BVDV)	DGHM/VAH 09/2001 (corr. EN 13727, EN 13624, EN 13897) EN 1040, EN 1275, EN 1278 DVV / RKI*	1 min.
Mycobacteria (Tb B)	DGHM 09/2001 (corr. EN 14348)	30 sec.
Fungi (<i>A. brasiliensis</i>)	EN 13624	15 min.
Non enveloped viruses	DW/RKI/carrier test	10 min.
Avian Influenza A	DVV / RKI*	15 sec.
HBV (DHBV)	DVV / RKI*	15 sec.
Rotavirus	DVV / RKI*	1 min.
Adenovirus	DVV / RKI*	1 min.
Polyomavirus	DVV / RKI*	5 min.
Poliovirus	surface carrier test	10 min.
NORO-virus (FCV)	DVV / RKI*	2 min.
NORO-virus (MNV)	EN 14478	1 min.

Expert reports are available upon request

Varenr.	Meliseptol® HBV tissue	Salgsenhed
18706	Meliseptol® ,HBV Tissue dispenser box med 100 servietter	20 stk./karton
18707	Meliseptol® ,HBV Tissue, refill 100 servietter	20stk/karton